



NanoDrop™ LOW VOLUME AUTOMATED PIPETTING

CAPABLE


Individual volume control with aspiration and dispensing from 100nL to 500uL

VERSATILE

Ability to handle wide range of reagents and solvents through the FEP, SS, and sapphire wetted parts

FLEXIBILITY



8 or 16 channel options for 96 down to 1536-wells as well as non-standard labware



www.bionexsolutions.com

©2024 BioNex Solutions Inc.

LIT-50004 Rev G Effective 2024-05



General Specifications

Dispense Channels	8 or 16 (NanoDrop Express needed for 16)
Aspiration	Up to 1000uL with <3uL dead volume
Software	Nanobuilder, RS232, API for integrations
Accuracy	<5% CVs at 1uL (with microsolenoids)